


Anti-CST3 Antibody (1L285)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 14 kDa.
Clone:	1L285

Applications

Verified Activity: 1. Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P01034
Synonyms:	Post-gamma-globulin;Cystatin-3;Cystatin-C;CST3;Neuroendocrine basic polypeptide;γ-trace; Post-γ-globulin;Gamma-trace

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481