

Anti-CST3 Antibody (3R866)



Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 14 kDa.
Clone:	3R866

Applications

Verified Activity:

1. Western Blot analysis of Cystatin C protein using antibody diluted at 1:1000.
2. Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.



Application:	IHC, WB
Recommended	WB: 1:1000-2000; IHC: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Uniprot ID: P01034
Synonyms: γ -trace; Neuroendocrine basic polypeptide; CST3; Gamma-trace; Cystatin-C; Cystatin-3; Post- γ -globulin; Post-gamma-globulin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481