

Anti-Desmin Antibody (2U295)



Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 53 kDa.
Clone:	2U295

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Mouse Skeletal Muscle Tissue using Desmin Mouse mAb diluted at 1:200.
2. Immunohistochemical analysis of paraffin-embedded Mouse Heart Tissue using Desmin Mouse mAb diluted at 1:200.



Application:	IF,IHC,WB
Recommended	WB: 1:200-1000; IHC: 1:100-200; IF: 1:200-400

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Uniprot ID: P17661

Synonyms: Desmin;DES

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481