

Anti-EBFP Polyclonal Antibody 2

Product Details

Conjugation: Unconjugated
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of EBFP recombinant protein, using TMAC-01217 diluted at 1:5,000.



Application: IF,IHC,IP,WB
Recommended WB: 1:2000-5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Synonyms: EBFP;Enhanced Blue Fluorescent Protein;BFP;Blue Fluorescent Protein

Research Background

Enhanced Blue Fluorescent Protein (EBFP) emits a strong blue fluorescence. EBFP can be used as protein "tags" to study the sub-cellular localization of proteins and/or their translocation upon stimulation or as markers for transfection in transient and stable expression systems.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481