

Anti-Phospho-E-Catenin (Ser655, Thr658) Polyclonal Antibody

Product Details

Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 100 kDa.
Purification:	Affinity purification using immunogen.

Applications

Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human E-Catenin around the phosphorylation site of Ser655/Thr658
Antigen Species:	Human
Uniprot ID:	P35221
Synonyms:	p-E-Catenin (S655, T658);p-E-Catenin (Ser655, Thr658);E-Catenin (p-S655, T658);E-Catenin (p-Ser655, Thr658)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481