


Anti-ECFP Polyclonal Antibody 2

Product Details

Conjugation: Unconjugated
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of ECFP Recombinant protein, using TMAC-01222 diluted at 1) 1:5,000
2) 1:10,000.



Application: WB
Recommended WB: 1:2000-5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Research Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481