

Anti-ECFP Polyclonal Antibody

Product Details

Ig Type: IgG
Conjugation: Unconjugated
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of ECFP recombinant protein, using TMAC-01223 diluted at 1:5,000.



Application: IP,WB
Recommended WB: 1:2000-5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Synonyms: Cyan fluorescent protein;Enhanced Cyan fluorescent protein

Research Background

Cyan Fluorescent Protein (CFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481