

Anti-Phospho-EGFR (Tyr998) Polyclonal Antibody

Product Details

Reactivity: Human, Mouse
Conjugation: Unconjugated
Molecular Weight: Theoretical: 180 kDa.
Purification: Affinity purification using immunogen.

Applications

Application: WB
Recommended WB: 1:1000-3000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human EGF Receptor around the phosphorylation site of Tyr998
Antigen Species: human
Uniprot ID: P00533
Synonyms: EGFR (p-Y998); p-EGFR (Tyr998); p-EGFR (Y998); EGFR (p-Tyr998)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481