


Anti-Phospho-ITK (Tyr512) Polyclonal Antibody 2

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Actual: 72 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western blot analysis of 293 using p-Emt (Y512) antibody.



Application: ELISA,IHC,WB
Recommended WB: 1:500-2000; IHC: 1:100-300;ELISA: 1:5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human Emt around the phosphorylation site of Y512
Antigen Species: human
Uniprot ID: Q08881
Synonyms: p-ITK (Tyr512);ITK (p-Tyr512);ITK (p-Y512);p-ITK (Y512)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481