

Anti-Phospho-EphA2 (Tyr772) Polyclonal Antibody

Product Details

Reactivity: Human
Conjugation: Unconjugated
Molecular Weight: Theoretical: 108 kDa.
Purification: Affinity purification using immunogen.

Applications

Application: WB
Recommended WB: 1:1000-3000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human EphA2 around the phosphorylation site of Tyr772
Antigen Species: human
Uniprot ID: P29317
Synonyms: p-EphA2 (Tyr772);EphA2 (p-Y772);p-EphA2 (Y772);EphA2 (p-Tyr772)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481