


Anti-Phospho-HER3/ERBB3 (Tyr1197) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 160-170 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western Blot analysis of 293T cells using Phospho-ErbB-3 (Y1197) Polyclonal Antibody.



Application: ELISA,WB
Recommended WB: 1:500-2000; ELISA: 1:40000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human ErbB-3 around the phosphorylation site of Y1197
Antigen Species: human
Uniprot ID: P21860
Synonyms: p-HER3/ERBB3 (Y1197);HER3/ERBB3 (p-Tyr1197);p-HER3/ERBB3 (Tyr1197);HER3/ERBB3 (p-Y1197)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481