

Anti-EYFP Polyclonal Antibody

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western blot analysis of EYFP recombinant protein, using TMAC-01407 diluted at 1:5,000.



Application:	IF,IHC,IP,WB
Recommended	WB: 1:3000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Synonyms: EYFP;Enhanced Yellow florescent protein;YFP;Yellow florescent protein

Research Background

Yellow Fluorescent Protein (YFP) is a genetic mutant of green fluorescent protein (GFP) originally derived from the jellyfish *Aequorea victoria*. Its excitation peak is 514nm and its emission peak is 527nm. YFP is a useful tool in cell and molecular biology thanks to its properties useful for fluorescence microscopy.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481