

Anti-Phospho-EZH2 (Thr487) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 85 kDa.
Purification:	Affinity purification using immunogen.

Applications

Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human EZH2 around the phosphorylation site of Thr487
Antigen Species:	Human
Uniprot ID:	Q15910
Synonyms:	p-EZH2 (Thr487);EZH2 (p-T487);EZH2 (p-Thr487);p-EZH2 (T487)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481