

Anti-Phospho-FGFR1/2/3/4 (Tyr653, 654) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 92 kDa.
Purification:	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Applications

Application:	WB
Recommended	WB: 1:1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human FGFR1/2/3/4 around the phosphorylation site of Tyr653+Tyr654
Antigen Species:	human
Uniprot ID:	P11362
Synonyms:	FGFR1/2/3/4 (p-Y653, 654);p-FGFR1/2/3/4 (Tyr653, 654);FGFR1/2/3/4 (p-Tyr653, 654);p-FGFR1/2/3/4 (Y653, 654)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481