


Anti-Fumarase Antibody (4Z791)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 50 kDa.
Clone:	4Z791
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of 1) 293T, 2) HepG2, 3) Hela, diluted at 1:3000.



Application:	IF,WB
Recommended	WB: 1:3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID:	P07954
Synonyms:	Fumarase;FH;Fumarate Hydratase Mitochondrial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481