

Anti-FLOT1 Antibody (60910)


Product Details

Reactivity:	Rat,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Actual: 49 kDa.
Clone:	60910
Purification:	Affinity purification using immunogen.

Applications

Verified Activity:

1. Western blot analysis of 1) Mouse Brain Tissue , 2) Rat Brain tissue, using TMAC-01514 diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Uniprot ID: O75955

Research Background

Two flotillin members have been described, flotillin-1 and flotillin-2. May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481