

## Anti-FMR1 Antibody (7X357)




### Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Clone:	7X357
Purification:	Affinity purification

### Applications

#### Verified Activity:

1. Western Blot analysis using FMR1 Monoclonal Antibody against Jurkat (1) and K562 (2) cell lysate.
2. Immunohistochemistry analysis of paraffin-embedded human brain tissues.
3. Immunofluorescence analysis of NIH/3T3 cells using FMR1 Monoclonal Antibody (green).



Application: ELISA,IF,IHC,WB

Recommended WB: 1:500-2000; IHC: 1:200-1000; IF: 1:200-1000; ELISA: 1:10000

**Properties**

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

---

**Antigen Details**

Immunogen: Recombinant Protein: human FMR1

Antigen Species: human

Gene ID: 2332

Uniprot ID: Q06787

Synonyms: Fragile X mental retardation protein 1;Protein FMR-1;FMR 1;FMRP

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481