


## Anti-Phospho-FOXO3 (Ser574) Polyclonal Antibody

### Product Details

Ig Type: IgG  
Reactivity: Human,Mouse,Rat  
Conjugation: Unconjugated  
Molecular Weight: Theoretical: 97 kDa.

### Applications

Verified Activity: 1. Western blot analysis FoxO3a (Phospho-Ser574) using TNF- $\alpha$  treated 293 whole cell lysates.



Application: WB  
Recommended WB: 1:1000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: A synthesized phosphopeptide: human FoxO3a around the phosphorylation site of Ser574  
Antigen Species: human  
Uniprot ID: O43524  
Synonyms: FOXO3 (p-Ser574);p-FOXO3 (S574);FOXO3 (p-S574);p-FOXO3 (Ser574)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481