

## Anti-LGALS3 Antibody (1V601)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 26 kDa.
Clone:	1V601

### Applications

Verified Activity:	Immunohistochemical analysis of paraffin-embedded human-colon using antibody diluted at 1:50.
Application:	IHC
Recommended	IHC: 1:100-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P17931
Synonyms:	L29;MAC2GAL3;AGE-R3;GALBPCBP35;GAL3;Galectin3;CBP35;Lectin L-29;Gal-3;LGALS3;Mac-2;L31;LGALS2;GALIG;galactin-3

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481