


Anti-LGALS3 Antibody (2I204)

Product Details

Ig Type: IgG1
Reactivity: Human,Mouse
Conjugation: Unconjugated
Molecular Weight: Theoretical: 26 kDa.
Clone: 2I204

Applications

Verified Activity: Western blot analysis of 1) MCF7, 2) 3T3 with Galectin-3 Mouse mAb diluted at 1:2,000.



Application: WB
Recommended WB: 1:2000-5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Uniprot ID: P17931
Synonyms: Mac-2;CBP35;galactin-3;MAC2GAL3;L29;GALIG;AGE-R3;Galectin3;GALBPCBP35;LGALS3;Gal-3;Lectin L-29;GAL3;L31;LGALS2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481