


Anti-GAPDH Antibody (4E643)

Product Details

Reactivity:	Human,Mouse,Rat,Monkey,Dog,Chicken,Hamouseter,Rabbit,Pig,Sheep
Conjugation:	Unconjugated
Molecular Weight:	Actual: 37 kDa.
Clone:	4E643
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of Hela, diluted at 1:3000.



Application:	WB
Recommended	WB: 1:2000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID:	P04406
Synonyms:	GAPD;glyceraldehyde-3-phosphate dehydrogenase;G3PD;HEL-S-162eP

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481