

Anti-GAPDH Antibody (9X189)




Product Details

Ig Type:	IgG1
Reactivity:	Human, Mouse, Rat, Monkey, Dog, Chicken, Hamouseter, Rabbit, Pig, Sheep, Insect, Yeast
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 37 kDa.
Clone:	9X189

Applications

Verified Activity:

- Western blot analysis of Hela (1), Rat brain (2), Rabbit Muscle (3), Sheep Muscle (4), and Mouse brain (5), diluted at 1:10000.
- The picture was kindly provided by our customer.
- The picture was kindly provided by our customer.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide

Uniprot ID: P04406

Synonyms: GAPD;HEL-S-162eP;G3PD;glyceraldehyde-3-phosphate dehydrogenase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481