



Anti-GFP Antibody (2Z569)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 2Z569

Applications

Verified Activity: 1. Western blot analysis of GFP transfected Hela, diluted at 1) 1:5000 2) 1:10000, 1) Input (control) 2) IP products, antibody dilution 1:200.



Application: IP, WB
Recommended WB: 1:500-2000; IP: 1:200

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Synonyms: GFP Protein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481