


## Anti-GFP Antibody (5E625)

### Product Details

Ig Type: IgG1  
Conjugation: Unconjugated  
Clone: 5E625

### Applications

Verified Activity: 1. 1) Input: HeLa Cell Lysate 2) IP product: IP dilute 1:200.



Application: IP, WB  
Recommended WB: 1:500-2000; IP: 1:200

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: Recombinant Protein  
Synonyms: GFP Protein

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481