


Anti-GFP Polyclonal Antibody 3

Product Details

Conjugation: Unconjugated
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of GFP Recombinant protein, using TMAC-01652 diluted at 1) 1:5,000, 2) 1:10,000.



Application: IF,IHC,IP,WB
Recommended WB: 1:5000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Synonyms: GFP Protein

Research Background

The green fluorescent protein (GFP) is a protein composed of 238 amino acid residues (26.9kDa) that exhibits bright green fluorescence when exposed to light in the blue to ultraviolet range. In cell and molecular biology, the GFP gene is frequently used as a reporter of expression.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481