


Anti-GFP Antibody (5I373)

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 26 kDa.
Clone:	5I373
Purification:	ProA affinity purified

Applications

Verified Activity: 1. Western blot analysis of GFP on K562 cells transfected with GFP using anti-GFP antibody at 1/100,000 dilution.



Application: ICC/IF,IHC,WB
Recommended WB: 1:10000-100000; IHC: 1:500-1000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Synonyms: GFP Protein

Research Background

The green fluorescent protein (GFP) was originally identified as a protein involved in the bioluminescence of the jellyfish *Aequorea victoria*. GFP cDNA produces a fluorescent product when expressed in prokaryotic cells, without the need for exogenous substrates or cofactors, making GFP a useful tool for monitoring gene expression and protein localization in vivo. Several GFP mutants have been developed, including EGFP, which fluoresce more intensely than the wildtype GFP and have shifted excitation maxima, making them useful for FACS and fluorescence microscopy as well as double-labeling applications. GFP is widely used in expression vectors as a fusion protein tag, allowing expression and monitoring of heterologous proteins fused to GFP.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481