


Anti-GFP-Tag Polyclonal Antibody

Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. 1ug GFP fusion protein+ Primary antibody dilution at 1) 1:5000 2) 1:20000. Secondary antibody was diluted at 1:20000.



Application:	IF, WB
Recommended	WB: 1:5000; IF: 1:50-200

Properties

Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: GFP-Tag
Synonyms:	GFP Protein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481