

## Anti-GSK3B Antibody (2C423)




### Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 46 kDa.
Clone:	2C423

### Applications

#### Verified Activity:

1. Western blot analysis of 1) Hela Cell Lysate, 2) 3T3 Cell Lysate, 3) Rat Brain Tissue Lysate using GSK3 $\beta$  Mouse mAb diluted at 1:1000.
2. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using GSK3 $\beta$  Mouse mAb diluted at 1:200.
3. Immunohistochemical analysis of paraffin-embedded Human Stomach Carcinoma Tissue using GSK3 $\beta$  Mouse mAb diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:1000-2000; IHC: 1:100-200

**Properties**

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

---

**Antigen Details**

Immunogen: A synthesized peptide

Uniprot ID: P49841

Synonyms: glycogen synthase kinase 3 beta;GSK-3 $\beta$ ;glycogen synthase kinase 3  $\beta$ ;GSK3 $\beta$

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481