



Anti-Phospho-GSK3B (Ser9) Polyclonal Antibody 2





Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 46 kDa.
Purification:	Antibodies were produced by immunizing rabbits with synthetic phosphopeptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific phosphopeptide. Non-phospho specific antibodies were removed by chromatography using non-phosphopeptide.

Applications

- Verified Activity:
1. Western blot analysis of extracts from 293 cells untreated (lane 1) or treated with serum (lane 2) using GSK3 β (Phospho-Ser9) Antibody TMAC-01736.
 2. Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using GSK3 β (Phospho-Ser9) Antibody TMAC-01736 (left) or the same antibody preincubated with blocking peptide (right).
 3. Immunofluorescence staining of methanol-fixed Hela cells showing cytoplasmic staining using GSK3 β (Phospho-Ser9) Antibody TMAC-01736.
 4. Immunohistochemical analysis of paraffin-embedded human Lung carcinoma tissue using GSK3 β (Phospho-Ser9) Antibody TMAC-01736.
 5. Immunofluorescence staining of methanol-fixed Hela cells showing cytoplasmic staining using GSK3 β (Phospho-Ser9) Antibody TMAC-01736.
 6. Western blot analysis of extracts from Hela cells, treated with EGF or calf intestinal phosphatase (CIP), using GSK3 β (Phospho-Ser9) Antibody TMAC-01736.





Application: IF,IHC,WB

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Synthesized phospho-peptide around the phosphorylation site of human GSK3β (phospho Ser9)

Antigen Species: Human

Uniprot ID: P49841

Synonyms: GSK3B (p-S9);p-GSK3B (S9);GSK3B (p-Ser9);p-GSK3B (Ser9)

Research Background

Participates in the Wnt signaling pathway. Implicated in the hormonal control of several regulatory proteins including glycogen synthase, MYB and the transcription factor JUN. Phosphorylates JUN at sites proximal to its DNA-binding domain, thereby reducing its affinity for DNA. Phosphorylates MUC1 in breast cancer cells, and decreases the interaction of MUC1 with CTNNB1/beta-catenin. Phosphorylates CTNNB1/beta-catenin.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481