


Anti-GST Antibody (8F169)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 8F169

Applications

Verified Activity: 1. Western blot analysis of 0.2ug GST protein, diluted at 1:5000.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Synonyms: GST

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481