

Anti-GST Tag Polyclonal Antibody 2


Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Molecular Weight:	Actual: 63 kDa.
Purification:	Antigen affinity purification.

Applications

Verified Activity:

1. Gel: 10%SDS-PAGE, Lysate: 2ug Recombinant GST tag protein, Primary antibody: 1/1000 dilution, Secondary antibody dilution: 1/8000, Exposure time: 15 seconds.



Application:	WB
Recommended	WB: 1:1000-8000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Full length fusion protein
Synonyms:	GST

Research Background

GST (Glutathione S-Transferase) is a 26kDa protein which was widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins. GST isoenzymes appear to play a central role in the parasite detoxification system. Other functions are also suspected including a role in increasing the solubility of haematin in the parasite gut.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481