


Anti-GST Tag Antibody (7C710)

Product Details

Conjugation: Unconjugated
Clone: 7C710
Purification: Affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of GST-tag fusion protein, using GST-Tag Mouse Monoclonal antibody (lane 1: 1:2000; lane 2: 1:5000).



Application: WB
Recommended WB: 1:3000-10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Full-length GST protein
Synonyms: GST

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481