

Anti-HA Tag Antibody (9C918)


Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Clone:	9C918
Purification:	ProA affinity purified

Applications

Verified Activity:

1. Western blot analysis of HA tag on HA tag recombinant protein lysates using anti-HA tag antibody at 1/1,000 dilution.



Application:	WB
Recommended	WB: 1:1000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Synonyms:	Hemagglutinin HA2;HA2;HA1;Hemagglutinin HA1 chain;HA epitope tag;hemagglutinin;HA-tag

Research Background

HA tag recombinant protein.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481