

Anti-HA Tag Antibody (9U839)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 9U839

Applications

Verified Activity: 1. Western blot analysis of HA Recombinant protein, diluted at 1:2000.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide
Synonyms: hemagglutinin; Hemagglutinin HA2; HA-tag; HA1; Hemagglutinin HA1 chain; HA2; HA epitope tag

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481