


Anti-HA Tag Antibody (7H489)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 7H489

Applications

Verified Activity: 1. Western blot analysis of Recombinant protein, diluted at 1) 1:5000 2) 1:10000.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide
Synonyms: HA1;hemagglutinin;HA-tag;Hemagglutinin HA2;HA epitope tag;Hemagglutinin HA1 chain;HA2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481