

Anti-Phospho-HDAC7 (Ser486) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 103 kDa.
Purification:	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.

Applications

Application:	IHC,WB
Recommended	WB: 1:500-2000; IHC: 1:50-200

Properties

Stability & Storage:	Store at -15°C to -25°C for 12 months (Do not lower than -25°C). Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human HDAC7A around the phosphorylation site of Ser486
Antigen Species:	Human
Uniprot ID:	Q8WUI4
Synonyms:	p-HDAC7 (S486);p-HDAC7 (Ser486);HDAC7 (p-Ser486);HDAC7 (p-S486)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481