


Anti-HER2/ERBB2 Antibody (6L396)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 138 kDa.
Clone:	6L396
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

- Verified Activity:
1. Western blot analysis of 1) Hela, 2) Mouse Brain, diluted at 1:4000.
 2. IHC staining of human breast cancer tissue, diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:2000-4000; IHC: 1:200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID: P04626

Synonyms: herstatin;HER-2;CD340;MLN 19;HERV7Q;HERVW;ERBB2;EGFR2;ENV;MLN19;ERVWE1;NGL;TKR1; ENVW;NEU;HER2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481