

Anti-ARNT Antibody (8R720)



Product Details

Ig Type:	IgG1
Reactivity:	Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 87 kDa.
Clone:	8R720

Applications

Verified Activity:

1. Western blot analysis of Mouse Brain Tissue with HIF-1 β /ARNT Mouse mAb diluted at 1:2,000.
2. Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using HIF-1 β /ARNT Mouse mAb diluted at 1:200.



Application:	IHC,WB
Recommended	WB: 1:1000-2000; IHC: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Uniprot ID: P27540

Synonyms: HIF-1-b;ARNT;HIF1-b

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481