


Anti-His Tag Antibody (7D784)

Product Details

Ig Type: IgG1
Conjugation: Unconjugated
Clone: 7D784

Applications

Verified Activity: 1. Western blot analysis of His-Tag recombinant protein, diluted at 1:5000 1) His-Tag transfected 293T 2) His-MBP recombinant protein.



Application: WB
Recommended WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481