

Anti-HDAC3 Antibody (4Q799)

Product Details

Reactivity: Human,Mouse,Rat
Conjugation: Unconjugated
Molecular Weight: Theoretical: 48 kDa.
Clone: 4Q799
Purification: Peptide affinity purified

Applications

Application: ICC,IHC,WB

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: Peptide
Uniprot ID: O15379
Synonyms: Histone deacetylase 3;HD3;Protein deacetylase HDAC3;HDAC3;RPD3-2;Protein deacetylase HDAC3;SMAP45

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481