


## Anti-Phospho-HIST1H1A (Thr17) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 17-25 kDa.
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

## Applications

Verified Activity: 1. Western blot analysis of extracts from Hela cells, untreated (-) or treated, 1:5000.



Application:	WB
Recommended	WB: 1:1000-2000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	Phospho peptide at Th17 of Histone H3
Uniprot ID:	Q02539
Synonyms:	HIST1H1A (p-Thr17);p-HIST1H1A (T17);HIST1H1A (p-T17);p-HIST1H1A (Thr17)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481