


Anti-Phospho-HIST1H1B (Thr17) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Monkey
Conjugation:	Unconjugated
Molecular Weight:	Actual: 31 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western blot analysis of SKOV3 using p-Histone H1 (T17) antibody.



Application: ELISA,IF,IHC,WB
Recommended WB: 1:500-2000; IHC: 1:100-300; IF: 1:200-1000; ELISA: 1:10000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human Histone H1 around the phosphorylation site of T17
Antigen Species: Human
Uniprot ID: P16401
Synonyms: HIST1H1B (p-T17);p-HIST1H1B (Thr17);HIST1H1B (p-Thr17);p-HIST1H1B (T17)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481