

## Anti-Histone H2B type 1-B Antibody (2K261)


### Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 14 kDa.
Clone:	2K261
Purification:	Affinity purification using immunogen.

### Applications

#### Verified Activity:

1. Western blot analysis of 1) HeLa, 2) 3T3, 3) Raw264.7, 4) Rat Brain, 5) Rat Kidney, using TMAC-01863 diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:1000-3000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

### Antigen Details

Uniprot ID: P33778

### Research Background

Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481