

Anti-Acetyl-Histone H2B type 1-B (Lys20) Polyclonal Antibody


Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 14 kDa.
Purification:	Affinity purification using immunogen.

Applications

Verified Activity:

1. Western blot analysis of extracts from Hela cells, untreated (-) or treated with TSA (1 μ M, 18 hr; +), using TMAC-01868 diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized acetylpeptide: human histone H2B around the acetylation site of Lys20
Antigen Species:	Human
Uniprot ID:	P33778
Synonyms:	Ac-H2BK20;Histone H2BK20-acetylated;Ac-Histone H2B type 1-B (K20);H2BK20ac;Ac-Histone H2B type 1-B (Lys20);Acetyl-Histone H2B type 1-B (K20)

Research Background

Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481