


## Anti-Tri-methyl-HIST1H3A (Lys36) Antibody (7X387)

## Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 15 kDa.
Clone:	7X387
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

## Applications

Verified Activity: 1. Western blot analysis of Hela, diluted at 1) 1:1000 2) 1:5000.



Application:	WB
Recommended	WB: 1:1000-3000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	tri methyl peptide at K36 of Histone H3
Uniprot ID:	P68431
Synonyms:	Tri-Me-HIST1H3A (K36);H3K36me3;TriMe-H3K36;Tri-methyl-HIST1H3A (K36);Tri-Me-HIST1H3A (Lys36);Histone H3K36-trimethylated

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481