


Anti-Tri-methyl-HIST1H3A (Lys4) Antibody (5G654)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 15 kDa.
Clone:	5G654
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of 1) HeLa, 2) Rat Testis tissue, 3) Raw264.7, diluted at 1:1000.



Application:	WB
Recommended	WB: 1:500-1000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	tri methyl peptide at K4 of Histone H3
Uniprot ID:	P68431
Synonyms:	H3K4me3;Tri-methyl-HIST1H3A (K4);Tri-Me-HIST1H3A (K4);Histone H3K4-trimethylated;TriMe-H3K4;Tri-Me-HIST1H3A (Lys4)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481