


Anti-Tri-methyl-HIST1H3A (Lys79) Antibody (9Q142)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 15 kDa.
Clone:	9Q142
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of Hela, diluted at 1:1000.



Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen:	tri methyl peptide at K79 of Histone H3
Uniprot ID:	P68431
Synonyms:	Tri-Me-HIST1H3A (K79);Histone H3K79-trimethylated;Tri-Me-HIST1H3A (Lys79);H3K79me3;Tri-methyl-HIST1H3A (K79);TriMe-H3K79

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481