


Anti-Tri-methyl-HIST1H3A (Lys9) Antibody (5M584)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 15 kDa.
Clone:	5M584
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of 1) HeLa, 2) Rat Testis tissue, 3) Raw264.7, diluted at 1:2000.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	tri methyl peptide at K9 of Histone H3
Uniprot ID:	P68431
Synonyms:	H3K9me3;Tri-Me-HIST1H3A (K9);Histone H3K9-trimethylated;TriMe-H3K9;Tri-Me-HIST1H3A (Lys9);Tri-methyl-HIST1H3A (K9)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481