


## Anti-Mono-methyl-HIST1H3A (Arg2) Polyclonal Antibody

### Product Details

Ig Type: IgG  
Reactivity: Mouse  
Conjugation: Unconjugated  
Molecular Weight: Theoretical: 15-17 kDa.

### Applications

Verified Activity: 1. Western blot analysis of 1) Raw264.7, diluted at 1:1000.



Application: WB  
Recommended WB: 1:500-1000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen: A synthesized monomethylpeptide: Histone H3 around the Monomethylation site of Arg2  
Synonyms: Mono-methyl-HIST1H3A (R2); MonoMe-H3R2; Histone H3R2-monomethylated; Mono-Me-HIST1H3A (Arg2); H3R2me; Mono-Me-HIST1H3A (R2)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481