


Anti-Mono-methyl-HIST1H3A (Lys79) Antibody (2J611)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 15 kDa.
Clone:	2J611
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity: 1. Western blot analysis of Hela, diluted at 1) 1:1000 2) 1:3000.



Application:	WB
Recommended	WB: 1:500-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	mono methyl peptide at K79 of Histone H3
Uniprot ID:	P68431
Synonyms:	Histone H3K79-monomethylated; Mono-Me-HIST1H3A (K79); H3K79me; Mono-Me-HIST1H3A (Lys79); Mono-methyl-HIST1H3A (K79); MonoMe-H3K79

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481